
VH SonicMist Ultrasonic Nebulizer System

Model 11-520


**Please read this instruction manual
completely before operating this unit.**

INSTRUCTION MANUAL • ENGLISH & ESPAÑOL


Stop!
Please ensure you have all of the following components before using your Nebulizer.


UltraSonic Unit


Mouthpiece #1


Nasal Piece


Child and Adult Masks

Mouthpiece #2

Extension Tube


Air Filter


5 Medication Cups

Unit Base

Battery Pack


Auto Adapter

AC Adapter

Storage Bag

Instruction Manual
FAQ Information Sheet
Replacement Parts
Order Form

IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.

INDEX

Care & Safety Information	2-3
Introduction.....	4
Product Features.....	7
Nebulizer Mask/Mouthpiece Selection.....	8
Power Option Use & Selection	9-10
Using Your Nebulizer	11-14
Cleaning & Disinfection	15-16
Changing the Filter	17
Troubleshooting	18
Warnings & Cautions	19
Device Symbols.....	20
Electromagnetic Compatibility	21-24
Product Specifications.....	25
Warranty	26
Instrucciones en Español.....	27-52

Toll-Free Customer Care Help Line: 1-866-326-1313
Monday – Friday 8:30 a.m. – 4:30 p.m. CST

Manufactured for:
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031
www.veridianhealthcare.com

Made in China

#93-1069 09/20
©2020 Veridian Healthcare, LLC

IMPORTANT SAFEGUARDS

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product. Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your nebulizer may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the nebulizer's efficacy and long life.

GENERAL CAUTIONS AND WARNINGS

- Adult supervision is recommended for use with children
- To avoid electrical shock keep unit away from water, never immerse the adapter or unit in liquid; do not reach for a unit that has fallen into water—immediately unplug the unit
- Do not handle the unit or power cord(s) with wet hands
- Do not use the unit if it has any damaged parts (including power cord or plug), if it has been submersed in water or dropped
- The unit should not be used where flammable gas, oxygen or aerosol spray products are being used
- If any abnormality occurs, discontinue use until the unit has been examined and repaired
- Do not leave the unit unattended when plugged in or operating
- Do not shake the unit when in operation
- Opening the cover during operation may cause injury
- Disconnect the unit from the electrical outlet before cleaning, filling and after each use
- Only use included attachments and accessories; do not use attachments from other nebulizer models
- Masks and mouthpieces should be replaced at least once a year
- Do not disassemble the unit or attempt to repair it

OPERATING CAUTIONS AND WARNINGS

- Do not allow nebulized medication to come in contact with eyes; nebulizer treatments are intended to be inhaled orally or nasally only

IMPORTANT SAFEGUARDS

CARE & SAFETY INFORMATION

- Change the water after each use
- The maximum capacity of the medicine cup is 7.5 ml.; do not overfill
- Always follow physician instruction for volume of medication used per treatment and necessary treatment time
- Never use this unit while operating a vehicle
- If you experience any discomfort or an abnormality occurs, stop using the unit immediately
- The unit should not be used for more than 30 minutes at a time, after use wait 30 minutes before beginning a new treatment

STORAGE CAUTIONS AND WARNINGS

- Never store the unit in direct sunlight, high temperature or humidity
- Keep the unit out of reach of small children
- Always keep the unit unplugged while not in use; it is also suggested that the battery pack be disconnected from the unit when not in use

CLEANINGS CAUTIONS AND WARNINGS

- Always check the selected mouthpiece, medication cup and air filter before each use; dirty or worn parts should be cleaned or replaced as necessary
- Never immerse the unit in water to clean as it may damage the unit
- Disconnect the unit from the electrical outlet or remove battery pack before cleaning
- Clean all necessary parts after each treatment as instructed in this guidebook; if the unit has not been used for an extended period, all parts should be cleaned prior to use

BATTERY CAUTIONS AND WARNINGS

- Included battery pack is suitable for use with VH SonicMist only, do not use with other units
- Do not use any other power source except those provided; unsupported power sources may damage unit
- Do not store the battery in direct sunlight, high temperature or high humidity
- Avoid placing the battery on metal or wet surfaces or inside a metal box; this may short-circuit the battery's positive and negative contact causing the battery to heat or explode

INTRODUCTION


Thank you for purchasing this Ultrasonic Nebulizer. Please read this instruction manual completely before using your unit to ensure safe and effective use, as well as long life of the product.

Ultrasonic nebulizers are used to efficiently deliver medication in the form of a mist to the lungs. Prescribed medicine is placed in the nebulizer and the ultrasonic transducer converts the liquid medication into an aerosol mist which is then inhaled. This method of delivery is beneficial in the treatment of asthma and other respiratory ailments, as the medicine is administered directly to the bronchial passages of the lungs to provide quick and effective relief. With proper use and care, your new nebulizer will provide you with many years of safe, effective breathing therapy treatments.

Federal law allows nebulizers to be sold by a physician's prescription only. These devices should be used only at the direction of a doctor or healthcare professional and should be used only as intended and described in this manual.

PRODUCT FEATURES

Illustrations of the nebulizer mechanics and detailed use instructions are included throughout this manual.

NEBULIZER MASKS/MOUTHPIECE SELECTION

This unit includes five styles of masks and mouthpieces, allowing for user preference in their desired medication delivery method. Replacement parts are available for order on the enclosed order form, by calling Customer Care at 866-326-1313 or at your local medical supply dealer.


Face mask (adult and child size included): face masks require use of the included extension tube. Select the appropriate size for the patient - child or adult - and place over the patients nose and mouth. Secure in place with the elastic strap around the patient's head and hold the nebulizer in place with a free hand.

Mouthpiece style #1: put the mouthpiece into the mouth, holding gently between the teeth and close the lips around the mouthpiece; hold the nebulizer in place with a free hand.

Mouthpiece style #2: this style requires use of the extension tube. Put the mouthpiece into the mouth, holding gently between the teeth and close the lips around the mouthpiece; hold the nebulizer in place with a free hand.

Nasal piece: gently insert the nasal pieces into the nose; hold the nebulizer in place with a free hand.

Proper use and assembly of each mouthpiece is further explained later in the 'USING YOUR NEBULIZER' section of this manual.


Adult Face Mask

Mouthpiece Style #1


Mouthpiece Style #2


Nasal Piece

POWER OPTIONS USE & SELECTION

This unit includes three convenient power options:

Rechargeable
Battery Pack

Auto Adapter

DC Auto
Adapter




Replacement power option parts are available for order on the enclosed order form, by calling Customer Care at 866-326-1313 or at your local medical supply dealer.

Rechargeable Battery Pack:

Note: Charge for 24 hours before first use; charge for 4-6 hours after initial charge. Battery can be charged with AC adapter plug or DC auto adapter.

1. Remove the battery pack from the housing.
2. Insert the adapter into the battery pack power port and plug the plug end into a wall socket.

A green light will illuminate on the battery and remain on while the unit is charging; this light does not necessarily indicate that the battery is fully charged, always follow the recommended charging times as noted above.

It is recommended that the battery be removed after each use and when the nebulizer will not be used for an extended period of time.

POWER OPTIONS USE & SELECTION

AC Adapter:

1. Be certain the battery pack is removed from the housing.
2. Insert the AC adapter directly into the underside of the housing, Fig 1. For hands-free operation, use with the included base, Fig 2.
3. Plug the AC adapter plug end into a wall socket.
4. Always unplug the adapter before cleaning and storage.


Fig. 1

DC Auto Adapter:

1. Be certain the battery pack is removed from the housing.
2. Insert the DC adapter directly into the underside of the housing, Fig 3.
3. Plug the adapter directly into the power source in the vehicle.
4. Always unplug the adapter before cleaning and storage.


Fig. 2


Fig. 3

Never leave the unit plugged in to a power source when not in use.

USING YOUR NEBULIZER

Before the first use of the nebulizer or if the unit has not been used for an extended period of time, it should be cleaned thoroughly as described in the **Cleaning and Disinfection** section of this manual.

The patient receiving treatment should be in a comfortable seated or slightly reclining position as treatments take several minutes.

1. Determine your power option and mouthpiece/mask style selection and gather your parts prior to starting the treatment.
2. Attach your selected power source.
3. Remove the clear nebulizer cover, Fig 1. Remove the medication cup if one is in place.
4. Fill the water reservoir with either tap or distilled water to the indicated level between the two lines within the reservoir, Fig 2. Use fresh water for every treatment.

NOTE: The recommended water temperature range is 68°F–104°F/20°C–40°C

The use of cool water (below 68°F/20°C) will result in a lower rate of nebulization.

The use of warm water (above 104°F/40°C) will result in a higher rate of nebulization.


Fig. 1


Fig. 2

USING YOUR NEBULIZER

5. Insert the medication cup as shown, Fig 1.
6. Pour the medication into the nebulizer cup following the physician prescribed dosage, Fig 2. The medication cup capacity is 7.5mL; overfilling the medication cup will reduce the nebulization rate.

Note: If your medication dosage volume is 2ml or less, or if the medication is very dense, the medication should be diluted with distilled water so that the volume is a minimum of 4 ml.

7. Replace the clear nebulizer cover, Fig 3.
8. Plug the power cord into the electrical outlet if necessary.

Fig. 1


Fig. 2



Fig. 3

Proper use of the mask or mouthpiece is important in helping to ensure effective and efficient medication delivery. Some pieces require an extension tube (included) for proper use. An extension tube is flexible and can be adjusted for a more comfortable position during treatment.

USING YOUR NEBULIZER

9. Attach the selected mouthpiece to the nebulizer and position the mouthpiece at the patient's mouth.

Remember that masks should cover the patient's mouth and nose; mouthpieces should be in the mouth with the lips closed around the end of the mouthpiece, Figs 4-7.


Face Mask

Fig. 4


Mouthpiece Style #1

Fig. 5

Mouthpiece Style #2

Fig. 6

Nasal Piece

Fig. 7

USING YOUR NEBULIZER

10. Press the **POWER** button on the nebulizer. Nebulization begins at a slower rate and will become stronger.
11. To increase or decrease the nebulization rate, press the **Fan** button until the desired level is reached. The lights on the housing will illuminate in correlation to the nebulization rate: 3 lights = highest rate, 1 light = lowest rate.
12. The nebulizer should remain vertical during treatment for proper nebulization of the medication and to prevent spilling of the liquid.
13. Ultrasonic treatment times average six minutes. The viscosity of the medication and airflow setting may result in longer treatment times.
14. When the treatment is complete, turn the unit off and unplug it from the outlet.
15. Always dispose of any remaining medication.
16. Allow the unit to cool. Disassemble, clean and store your nebulizer. Be certain that the nebulizer pieces are completely dry before storage.

It is important to complete the entire prescribed dose of medication for maximum effectiveness. Consult with your prescribing physician or healthcare professional for information on your specific medication therapy treatment plan.

CLEANING & DISINFECTION

The nebulizer and mouthpiece should be thoroughly cleaned after each use unless a different cleaning regimen is specified by your healthcare professional; additional hygienic precautions may be required for use with high-risk patients. Never attempt a subsequent medication treatment without first cleaning your unit and mouthpiece.

Cleaning Procedures:

1. Ensure that the power cord is disconnected from the power outlet and that the battery pack is disconnected from the nebulizer housing body.
2. Remove the clear lid from the nebulizer body, Fig 1.
3. Wash the mouthpiece, fan, medicine cup and air filter cover in warm water with mild detergent. Rinse thoroughly, Fig 2.
4. Allow all pieces to air dry before use or storage. Do not store any parts wet as remaining moisture creates a risk of bacterial growth.
5. Rinse the water reservoir with warm water, do not use detergent. As needed, wipe out the reservoir with a damp cloth. Do not replace the lid or medication cup until the unit is completely dry.
6. The body of the unit and the outside of the battery pack can be wiped down with a soft cloth; if necessary, use a damp soft cloth with a non-abrasive cleaner. Never immerse the nebulizer body or battery in liquid.


Fig. 1


Fig. 2

CLEANING & DISINFECTION

Disinfecting Procedures:

It is recommended that the parts be disinfected monthly if the unit is used on a regular medication therapy schedule.

1. Remove the clear lid from the nebulizer body.
2. Rinse the mouthpiece, fan, medicine cup and air filter cover in hot water. Mix one part white vinegar with three parts distilled water, in a volume sufficient to fully immerse the pieces.
3. Submerge pieces in the vinegar mixture for thirty minutes.
4. Rinse everything thoroughly with water.
5. Allow all pieces to air dry completely, or dry them thoroughly with a clean, soft towel. Do not store any parts wet as remaining moisture creates a risk of bacterial growth.

Never use soap pads or powder cleanser to clean any part of the unit. Never use any cleaning products that may be harmful when swallowed or inhaled.


It is recommended that mouthpieces and masks be replaced periodically depending on the frequency of therapy treatments or when damaged or cracked. To order replacement parts, use the included order form, call Customer Care at 866-326-1313 or contact your local dealer.

CHANGING THE FILTER

- Never operate the unit without a filter.
- Replace the filter after every ten uses, or sooner if it turns gray, appears contaminated or is blocked or becomes wet.
- Use only VH SonicMist filters; additional filters are available for order using the included order form, calling Customer Care at 866-326-1313 or contacting your local dealer.
- Do not substitute cotton or any other materials for the filters.
- Do not attempt to clean or reuse a filter.

To replace the filter, follow these steps:

1. Change the filter only when the unit is off, unplugged and completely cool.
2. Gently pull the filter cap out of the cover; if necessary, use a non-metal tool to help lift-up cap.
3. Remove the filter from the filter cap.
4. Place a new filter in the filter cap.
5. Return the filter cap to its place, ensuring that it is firmly in place.

TROUBLESHOOTING

The nebulizer produces no aerosol, or a small amount of aerosol.

- There may be too much or too little medicine in the nebulizer. The amount of liquid should be between 0.5 ml and 7.5 ml.
- The medication cup may be damaged; check and replace.
- The water reservoir may be dirty; wipe clean and attempt again.
- The nebulizer should be held vertically during use.
- The medication may not be suitable for nebulization. Check with your physician.
- The filter may be in need of replacement

Some medicine remains in the nebulizer at the end of my treatment.

This is normal. Never keep leftover medication.

Treatments with the mask takes longer than treatments with the mouthpiece.

This is normal. Less medication is delivered per breath with the mask because the aerosol is combined with ambient air through the holes in the mask.

Unit does not operate when the POWER button is pressed.

If using the AC adapter or auto adapter, check connections between the adapter and the unit and the adapter and the outlet. If using the battery, check that the battery is properly connected. Charge the battery.

Fan is not working properly.

Check that the fan is firmly pushed onto the motor shaft.

Fan speed light changes from red to green during use.

The battery needs to be charged. Nebulization will slow and then cease.

WARNINGS & CAUTIONS

Nebulizer therapy is intended for use with physician prescribed medication only.

Over the counter medication, essential/aromatherapy oils or alcohol are not suitable for inhaling using a nebulizer. Highly viscous substances may damage the unit and make it unsuitable for use. Never use nebulizer therapy for anything other than physician prescribed respiratory care therapy.

Use common sense electrical caution when operating the device. This nebulizer runs on DC power. Always turn off the power and unplug the device before cleaning or if the unit is not working properly.

DEVICE SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning
	Read This Manual—All included manuals should be read prior to first use.
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator.
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out it's intended function.
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal .
	Keep Dry—This device should be kept dry; never submerge the unit. Consult with the Care and Maintenance section of this manual for information on cleaning your device.
	Double Insulation—This unit is double-protected and thus corresponds to protection class 2.

ELECTROMAGNETIC COMPATIBILITY

Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emissions

This device is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Emission Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	This monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

ELECTROMAGNETIC COMPATIBILITY

Table 2

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic immunity

This device is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this monitor requires continued operation during power mains interruptions, it is recommended that this monitor be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

ELECTROMAGNETIC COMPATIBILITY

Table 3

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This device is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range.** Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this monitor.

** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

ELECTROMAGNETIC COMPATIBILITY

Table 4

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Recommended separation distances between portable and mobile RF communications equipment and this monitor

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PRODUCT SPECIFICATIONS

Name	VH SonicMist Ultrasonic Nebulizer
Model Number	11-520
Medication Capacity	Maximum 7.5 ml
Particle Size	0.5 to 10 μ m
Ultrasonic Frequency	2.5MHz+200kHz (-100 kHz)
Average Nebulization Rate	\geq 0.2 ml / minute
Operating Conditions	Temperature: 41oF – 104oF (5oC – 40oC); Humidity: 10% – 90%
Unit Dimensions (W x H x D)	3.35" x 2.35" x 5" (without battery) 3.35" x 2.35" x 6.25" (with battery)
Unit Weight	8 oz. (without battery) 13.8 oz. (with battery)
Power Connection	12V DC
Included Accessories	Rechargeable battery pack, 12-volt AC adapter, DC auto adapter, child mask, adult mask, extension tube, nasal piece, 2 mouthpiece styles, 5 medication cups, 1 air filter, base

Electronic compatibility in accordance with IEC 60601-1-2

Safety requirements in accordance with IEC 60601-1

Specifications are subject to change without notice.

3-YEAR LIMITED WARRANTY

The warrantor guarantees that its Ultrasonic Nebulizer will be free from manufacturing defects under normal use for a period of three years. This warranty covers only normal use, and does not apply to use in any clinical or commercial applications. This warranty does not cover batteries or other power sources that may be provided with or used with the nebulizer unit. This warranty is voided if the nebulizer unit is misused or abused in any manner.

If the nebulizer fails to operate within the warranty period while under the original purchaser owns the product, return it postage prepaid a to: **Veridian Healthcare, Attn: Repair Department, 1175 Lakeside, Gurnee, IL, 60031**; include your name, address, phone number and brief description of the problem with your nebulizer. Take care to carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. The warrantor will repair or replace the defective unit at our option. Repair or replacement of the defective unit, at the warrantor's option, is the sole remedy under this warranty.

ANY IMPLIED WARRANTIES WHICH THE PURCHASER MAY HAVE ARE LIMITED IN DURATION TO THE TIME THAT THE ORIGINAL CONSUMER PURCHASER OWNS THE PRODUCT. Some states do not allow limitations or how long an implied warranty lasts, so the above limitation may not apply to you.

This warranty constitutes the warrantor's only reasonability and obligation to repair and/or replace materials or components, or refund the purchase price. The warrantor will not be responsible for any indirect, incidental, special consequential or punitive damages or other loss, including, but not limited to damage to or loss of other property or liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchase price of the unit. Some states do not allow the exclusive or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

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