
Steam Inhaler and Beauty Mask


Model 11-525


**Please read this instruction manual
completely before operating this unit.**

INSTRUCTION MANUAL • ENGLISH & ESPAÑOL


STOP!
PLEASE ENSURE YOU HAVE ALL OF THE FOLLOWING
COMPONENTS BEFORE USING YOUR STEAM INHALER.


Steam Inhaler

Small Face Mask +
Large Facial Mask

Measuring Cup

Secondary Fill Tank

Instruction Manual

**IF YOU ARE MISSING ANY PARTS, INCLUDING
INSERTS OR INSTRUCTION MANUALS, DO
NOT RETURN TO PLACE OF PURCHASE.
CONTACT CUSTOMER CARE AT 866-326-1313.**

INTRODUCTION

Thank you for purchasing this Veridian Steam Inhaler. **PLEASE READ THIS INSTRUCTION MANUAL COMPLETELY BEFORE USING YOUR UNIT TO ENSURE SAFE AND EFFECTIVE USE, AS WELL AS LONG LIFE OF THE PRODUCT.**

The included two mask sizes allow this product to serve the dual purposes of hot steam inhalation therapy and facial beauty treatments. Steam inhalation therapy is a natural means for helping to relieve the congestion of bronchitis, the common cold, allergies and a variety of other respiratory and sinus ailments. Steam inhalation opens up congested sinuses and lung passages, allowing you to discharge mucus, breathe more easily and heal faster. Steam facials open the pores of the skin, helping to naturally remove dirt and residue build-up without the use of harsh soaps and chemicals.

INDEX

Introduction	2
Important Safeguards - Care & Safety Information	4-5
Product Features	6
Face Mask Selection	7
Using Your Steam Inhaler – Hot Steam Inhalation Therapy	8-10
Using Your Steam Inhaler – Facial	11
Cleaning and Disinfection	12-13
Troubleshooting	14
Electromagnetic Compatibility	15-18
Product Specifications	19
Warranty	20
Instrucciones en Español.....	21-40

Toll-Free Customer Care Help Line: 1-866-326-1313
Monday – Friday 8:30 a.m. – 4:30 p.m. CST

Manufactured for:
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031
www.veridianhealthcare.com

Made in China

#93-1070 05/23
©2017 Veridian Healthcare, LLC

IMPORTANT SAFEGUARDS

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product. **Failure to read and observe all precautions could result in personal injury or equipment damage.**

Steam output may reach in excess of 109°F (42.8°C). Always use caution when handling the unit. Not intended for use by children or infirmed persons without close adult supervision.

Improper care or use of your steam inhaler may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the steam inhaler's efficacy and long life.

GENERAL CAUTIONS AND WARNINGS

- Keep out of reach of children
- To avoid electrical shock keep the steam inhaler away from water, never immerse the unit in liquid; do not reach for a unit that has fallen into water—immediately unplug the unit
- Do not handle the unit or power cord with wet hands
- Do not use the unit if it has any damaged parts (including power cord or plug), if it has been submersed in water or dropped
- If any abnormality occurs, discontinue use until the unit has been examined and repaired
- Do not leave the unit unattended when plugged in or operating
- Do not shake the unit when in operation
- Opening the cover during operation may cause injury
- Disconnect the unit from the electrical outlet before cleaning, filling and after each use
- Only use included attachments and parts; do not use attachments from other brands or models
- Do not disassemble the unit or attempt to repair it

IMPORTANT SAFEGUARDS

CARE & SAFETY INFORMATION

OPERATING CAUTIONS AND WARNINGS

- Do not put masks up to eyes or allow steam to come in contact with eyes
- Change the water after each use
- Do not overfill the heating tank or secondary fill tank
- Only use clean pure water in the heating tank
- Never use this unit while operating a vehicle or in the bathtub or shower
- If you experience any discomfort or an abnormality occurs, stop using immediately
- The unit should not be used for more than 10 minutes at a time; always wait until the unit has completely cooled before beginning a treatment
- Do not block the side vents while in use
- Always turn off and unplug immediately after use


STORAGE CAUTIONS AND WARNINGS

- Never store the unit in direct sunlight, high temperature or humidity
- Keep the unit out of reach of small children
- Always keep the unit unplugged while not in use

CLEANING CAUTIONS AND WARNINGS

- Never immerse the unit in water to clean as it may damage the unit
- Disconnect the unit from the electrical outlet before cleaning
- Clean all necessary parts after each treatment as instructed in this instruction manual; if the unit has not been used for an extended period, all parts should be cleaned prior to use

PRODUCT FEATURES

DETAILS OF INTERIOR PARTS ARE INCLUDED THROUGHOUT THIS MANUAL.

FACE MASK SELECTION

This steam inhaler includes two mask styles—a small nose/mouth style and a large facial style. Detailed instructions for use of the different mask styles are found on pages 8 thru 11 of this manual.

The **small nose/mouth style mask** is suitable for use with hot steam inhalation therapy. The small cone of this mask helps direct the steam into the mouth and nose. The user's mouth and nose should fit comfortably into the mask; the notched portion at the top of the mask allows room for the bridge of the nose.

The **large facial style mask** is suitable for use with hot steam facials. The large cone of this mask allows the steam to penetrate the skin pores around the cheeks, mouth, nose and chin. The user's face EXCLUDING the eyes should fit comfortably into the mask. Never cover your eyes with the mask, even when keeping the eyes closed.

NEVER REST YOUR FACE ON THE MASK. THE MASKS ARE PLASTIC AND WILL BECOME PLIABLE WHEN HEATED BY THE STEAM; SUPPORTING THE HEAD/FACE ON THE MASK COULD CAUSE THE UNIT TO TIP OR FALL.

USING YOUR STEAM INHALER – HOT STEAM INHALATION THERAPY

Before proceeding, the individual receiving therapy should be seated in a comfortable position as treatments take several minutes. Keep in mind the height of the user and of the product to ensure the mask is at the desired level; the steam inhaler should never be lifted or handled while operating.

STEAM THERAPY TREATMENTS

1. The unit should be situated on a flat, stable surface within operating range of the user and an electrical outlet. Ensure that the air vents on the side of the unit housing are not blocked or obstructed.
2. Remove the top cover of the housing by pressing the ridged area on the underside of the extension tube and pulling up on the top cover, Fig. 1.
3. Open the water heater tank by turning the cover counter-clockwise, Fig. 2.
4. Using the included measuring cup, add 20 ml of pure tap or distilled water to the water heater tank, Fig. 3. Note the fill line in the tank. Do not overfill. Never add saline, essential oils or other additives as this will damage the unit.
5. Replace the water tank cover and close securely by turning clockwise.


Fig. 1


Fig. 2


Fig. 3

USING YOUR STEAM INHALER – HOT STEAM INHALATION THERAPY

6. Remove the secondary fill tank, Fig. 4.
7. Add 80 ml of tap/distilled water to the larger capacity side of the secondary fill tank (the side where the suction tube feeds), Fig. 5.
8. Replace the secondary fill tank, making sure the suction tube feeds into the larger side of the fill tank. The left/smaller side of the secondary fill tank is an overflow tank to collect water. Never place liquid in this side of the tank.
9. Check to make sure the nozzle and suction tube are in place and that the suction tube is not bent, Fig. 6. Reattach the top cover to the body of the unit.
10. Place the smaller, mouth/nose mask onto the extension tube and adjust the extension tube to the desired height, Fig. 7. The height/length of the tubing will not affect the steam output or treatment time.


Fig. 4


Fig. 5


Fig. 6


Fig. 7

USING YOUR STEAM INHALER – HOT STEAM INHALATION THERAPY

11. Make sure the dial switch on the power cord is in the 'OFF' position and plug the power cord into an electrical outlet. Turn the dial switch to the ON position, Fig. 1. The light on the front of the unit under the secondary fill tank will illuminate.
12. Bring face to the mask, Fig 2. Do NOT attempt to lift the unit to face. The face may safely make contact with the mask, but should not rest on or pull on the mask.
13. Unit will take approximately 60 seconds to boil water and begin producing steam.
14. Inhale and exhale normally. The particle size can be adjusted by sliding the particle adjustment slide switch, Fig. 3.

Up—produces larger particles, ideal for inhalation through the mouth

Down—produces smaller particles, ideal for inhalation through the nose

15. Continue treatment for as long as comfortable or until the steam lessens and quits, approximately 6-9 minutes.
16. When treatment is complete, turn off unit immediately at the dial switch on the cord and unplug the unit. Do not attempt to move or disassemble the unit until it has cooled completely, approximately 10 minutes.
17. Disassemble and clean your steam inhaler, Fig. 4. All parts must be dry before storage.


Fig. 1

Fig. 2


Fig. 3



Fig. 4

USING YOUR STEAM INHALER – FACIAL TREATMENTS

Steam facials help to open the pores in the skin and coax out the dirt, oil, makeup residue and impurities that have built up within the skin.

It is recommended that you seek the advice of a licensed esthetician or cosmetology professional for direction on the recommended frequency and duration of steam facials.

Use of the Steam Inhaler for a facial requires the same steps as when using this product for hot steam inhalation therapy, utilizing the larger facial mask rather than the small mask.

STEAM FACIAL TREATMENTS

1. Review steps 1-13 in the 'USING YOUR STEAM INHALER - HOT STEAM INHALATION THERAPY' section of this manual for guidelines on the set-up and operation of the unit.
2. Place the large facial mask on the unit and proceed with your treatment.
3. Inhale and exhale normally. Continue treatment for as long as comfortable or until the steam lessens and quits, approximately 6-9 minutes.
4. When treatment is complete, turn off unit immediately at the dial switch on the cord and unplug the unit. Do not attempt to move or disassemble the unit until it has cooled completely, approximately 10 minutes.
5. Disassemble and clean your steam inhaler. Be certain that all parts are dry before storage.


CLEANING AND DISINFECTION

The steam inhaler mask(s) and accessories should be thoroughly cleaned after each use and disinfected each day if being used to treat an illness unless a different cleaning regimen is specified by your healthcare professional.

Always make certain that the unit has cooled completely after a therapy session before attempting to dismantle or clean the steam inhaler.

CLEANING/RINSING:

1. Ensure that the power is set to OFF and the power cord is disconnected from the power outlet.
2. Disconnect the mask(s). Remove the secondary fill tank, the top portion of the unit housing and the lid of the water heating tank.
3. Empty any remaining water from the water heat tank.
4. Use warm water to rinse the mask(s) and secondary fill tank.
5. Allow the unit to completely air dry or dry thoroughly with a clean, soft towel. Do not store any parts wet as remaining moisture creates a risk for bacterial growth.
6. Never wash or submerge steam inhaler housing or power cord. Occasionally wipe the body of the unit with a damp cloth and allow it to air dry. Only clean the body of the unit when it is unplugged and completely cool.

NOTE: Steam inhaler parts and masks are NOT dishwasher safe.

DISINFECTION:

The same general process may be followed when disinfecting the parts of the steam inhaler. Mix one part white vinegar with three parts distilled water in a volume sufficient to fully immerse all of the parts. Submerge the parts in the water and vinegar mixture for thirty minutes, then rinse all parts thoroughly with water.

CLEANING AND DISINFECTION

CLEANING THE STEAM NOZZLE

Water with a high mineral count may clog the steam nozzle. If the steam inhaler begins to produce a significant less amount of steam or no steam but has power and water in the water heater tank, the steam nozzle may need to be cleaned. Be certain that the unit has entirely cooled before attempting to clean the steam nozzle.

1. Locate the small spring pin on the underside of the unit, Fig. 1. Use a non-metal tool to pop the pin out.
2. Remove the top housing of the unit. Locate the steam nozzle, twist counter-clockwise and lift up the nozzle as shown, Fig. 2.
3. Use the small spring to thread each of the two small holes in the steam nozzle to clear the holes, Fig. 3.
4. Replace the steam nozzle. Make sure the suction tubing is straight before attempting to use, Fig. 4.


Fig. 1


Fig. 2


Fig. 3


Fig. 4

TROUBLESHOOTING

POWER LAMP/LIGHT IS NOT ILLUMINATED

- Check that the power cord is connected to a working outlet
- Check that the power switch is in the ON position

LITTLE OR NO STEAM OUTPUT

- Check that the power is connected and switch is in the ON position
- Check the water level in the water heating and secondary fill tank and fill as necessary
- Clean the steam nozzle

WATER IS LEAKING

- Check that the secondary tank is not overfilled and is at the proper fill level
- Make sure the unit is on a flat surface and not tilted

STEAM IS LEAKING

- Ensure the water heating tank lid is firmly closed
- Ensure the top of the steam nozzle is firmly locked in place

ELECTROMAGNETIC COMPATIBILITY

Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an

Emission Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	This monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

ELECTROMAGNETIC COMPATIBILITY

Table 2

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this monitor requires continued operation during power mains interruptions, it is recommended that this monitor be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

ELECTROMAGNETIC COMPATIBILITY

Table 3

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range.** Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this monitor.

** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

ELECTROMAGNETIC COMPATIBILITY

Table 4

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Recommended separation distances between portable and mobile RF communications equipment and this monitor

This monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PRODUCT SPECIFICATIONS

NAME	Steam Inhaler and Beauty Mask
MODEL NUMBER	11-525
POWER SOURCE	120V/60Hz
POWER CONSUMPTION	150W
WATER HEATING TANK	25 ml capacity
SECONDARY FILL TANK	80 ml capacity
MEASURING CUP	25 ml capacity
TREATMENT/OPERATION TIME	6-9 minutes
UNIT DIMENSIONS	5-1/8" x 3-3/4" x 10"
UNIT WEIGHT	1 lb. 5 oz.
ACCESSORIES	Small mouth/nose style mask, large facial style mask, secondary fill tank, measuring cup

Specifications are subject to change without notice.

LISTED
Steam Inhaler
96GG
E239062

SAVE THESE INSTRUCTIONS

1-YEAR LIMITED WARRANTY

Veridian Healthcare guarantees that its Steam Inhaler will be free from manufacturing defects under normal use for a period of one-year. This warranty covers only normal use, and does not apply to use in any clinical or commercial applications. This warranty is voided if the unit is misused or abused in any manner.

If the steam inhaler fails to operate within the warranty period while under the original purchaser owns the product, return it postage prepaid to: **Veridian Healthcare, Attn: Repair Department, 1175 Lakeside, Gurnee, IL, 60031**; include your name, address, phone number and brief description of the problem with your nebulizer. Take care to carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. Veridian Healthcare will repair or replace the defective unit at our option. Repair or replacement of the defective unit, at the warrantor's option, is the sole remedy under this warranty.

ANY IMPLIED WARRANTIES WHICH THE PURCHASER MAY HAVE ARE LIMITED IN DURATION TO THE TIME THAT THE ORIGINAL CONSUMER PURCHASER OWNS THE PRODUCT. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

This warranty constitutes the warrantor's only responsibility and obligation to repair and/or replace materials or components, or refund the purchase price. Veridian will not be responsible for any indirect, incidental, special consequential or punitive damages or other loss, including, but not limited to damage to or loss of other property or liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchase price of the unit. Some states do not allow the exclusive or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

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